

## PRODUCT SPECIFICATIONS



- ✓ Covers the wrist gently and smoothly
- ✓ Can be used both hand and leg wrist
- ✓ With connection ribbons of 40 cm fix the wrist safely
- ✓ Absorbs sweat with its cotton fabric
- ✓ In case of long term use does not cause any allergy and irritation

- ✓ In case of long term use does not cause any allergy and irritation
- ✓ One fixation band in each pouch
- ✓ Band made of a master material, consists of 3 different material laminated to each other with the thickness of 4mm (+-0.5mm)
- ✓ Skin touch surface of the master material made of cotton fabric, upper part to cover the fixation velcro easily made of feathered jersey, middle part made of foam to provide elasticity and smoothness