

### PRODUCT SPECIFICATIONS



- ✓ The Crisis Set contains 16 pcs Steel Lock and 2 pcs Magnetic Key
- ✓ The set is placed into a protective bag
- ✓ The products are made of 100% terylene and polyester material
- ✓ Washable at 95 °C
- ✓ 2 years warranty

- ✓ Safebelt patient safety belts have a wide range to respond all cases of prevent and limited of undesired movements
- ✓ By providing a secure fixation offers a maximum safety both for patients and around
- ✓ Provides maximum freedom and comfort sensation to the patients where the patient movement is restricted or controlled



- ✓ Application is easy and practical
- ✓ It is used to prevent patient movement at the time of crisis
- ✓ The Crisis Set consists of;  
2 pcs Wrist Restraints Belt,  
2 pcs Ankle Restraints Belt, Waist Roll Belt,  
Shoulder Restraints,  
Pelvic Support Belt and Multi function extender belt