






PRODUCT SPECIFICATIONS



- Both sides micro laminated, unfeathered, unfuzzy, smooth and absorbative
- ✓ Hydrofiber nonwoven structured
- ✓ Product has got 3 sizes; small, medium and large

Elastic hole fits all models of tracheal cannulas.

Hydrophilic PU FOAM Dressing

- ✓ The skin-coloured cover of the dressing is semi-permeable and free of latex and adhesives
- ✓ Provides ideal removal of wound secretions as well as good cushioning
- → Absorbation capacity 30 times of its weight

