

PRODUCT SPECIFICATIONS



- ✓ Single use
- ✓ One band in each pouch
- ✓ Prevents from the UV rays up to 450 Nm
- ✓ Both type has 3 sizes; small, medium and large

Deluxe Model

- ✓ Y Shaped, it has a special design combination of the elastic, smooth band and UV protected eye band
- ✓ Supporter band covers the head safely and with its light permeability does not prevent the UV rays to the head
- ✓ Its protective part consists of carbon foam which protects the eyes %100 from the UV rays
- ✓ Air permeability texture of the product prevents the sweat

Classic Model

- ✓ Designed to protect the eye most ideal form
- ✓ Covers the face closely and softly
- ✓ Easy to use and fits the baby's head
- ✓ Extremely safe with connection ribbons
- ✓ Absorbs the sweat with its cotton fabric in case of long term use does not cause any irritation

